



Easy to transport and store, CB balance delivers performance with compact portability. Battery operation means it can be used anywhere, while auto power-off helps conserve power. A sealed keypad and one-piece housing are simple to clean. The grade 304 stainless steel pan can be wiped down quickly and accommodates either single items or bulk powders. A backlit display shows the user the weight and remaining battery level. An optional power adapter is available

## Convenient features for portable precision weighing



### Features

- Vivid, backlit LCD easily visible in any lighting conditions
- Large, grade 304 stainless steel pan allows swift cleaning
- Lightweight design simplifies transport and portability
- Sealed keypad protects against dirt and spills
- Level indicator and adjustable feet ensure proper balance setup for optimum weighing results
- Auto power-off to save energy
- External calibration allows for verification and adjustment with weights
- Simple operation allows quick startup and use with little or no training
- Zero-tracking feature ensures display returns to zero reading
- Low-battery indicator alerts users when batteries need replacement

### Applications

- Weighing

### About ADAM:

For more than 40 years, professionals worldwide have relied on Adam Equipment for an extensive selection of dependable and affordable weighing equipment. Adam provides the right balance of speed, performance and value.

Model	CB 501	CB 1001	CB 3000
Capacity	500g	1000g	3000g
Readability	0.1g	0.1g	1g
Pan Size	5.1" ø / 130mm ø		
Weighing Units	g, lb, oz		
Stabilization Time (sec)	2		
Calibration	External calibration		
Display	LCD with 0.5" / 12mm-high digits		
Power Supply	2 x AA batteries / Optional 6VDC 200mA adapter		
Operating Temperature	0° to 40°C / 32° to 104°F		
Housing	ABS plastic		
Overall Dim.	6.3"x8.1"x1.8" / 160x205x45mm (wxdxh)		
Net Weight	0.9lb / 0.4kg		

## Accessories

Item number	Description
700400120.....	6VDC 200mA adapter
700230015.....	In-use wet cover
700660290.....	Calibration certificate
700100060.....	M1 500g Calibration Weight
700100029.....	M1 1kg Calibration Weight

Supplied by